

# TIME LINE



We respectfully ask for the following time to complete your requested work.

Fabricate soft tissue models with analogs	2 Days
Bite rim	2 Days
Verification jig	2 Days
Articulation & set denture teeth for try-in	10 Days
Revise wax try-in	5 Days
Scan, mill and finish full arch zirconia	15 Days
Scan, mill and finish PMMA provisional	10 Days

All times are working days in our laboratory. This does not include the day received due to most deliveries arriving mid-afternoon. Times are exclusive of shipping/travel time.

